

RHEXIS

**STERILE SINGLE USE
INSTRUMENTATION
FOR OCULAR SURGERY**

Catalogue



RHEXIS CATARACT SURGERY

STERILE SINGLE USE INSTRUMENTATION FOR OCULAR SURGERY



SAFETY — **EFFECTIVENESS**
ECONOMY

TABLE OF CONTENTS CATARACT SURGERY

DESIGNATION	PAGE
● ANESTHESIA NEEDLES	4
● CYSTOTOMES	5
● HYDRODISSECTION CANNULAS	6
● IRRIGATION CANNULAS	7
● VISCO-ELASTIC CANNULA	8
● CORTEX ASPIRATION	9-11
● LACRIMAL CANNULA	12
● EYE SHIELD	13
● EYES SPONGES	14
● ABSORBENT STICK	15
● PRODUCT INFORMATION	16
● INDEX	17
● PRICE LIST	19



Designed and manufactured by : **LPH PHARMACEUTICAL**



PRODUCT INFORMATION

ANESTHESIA NEEDLES

Peribulbar, retrobulbar and peri-retrobulbar needles, ready for use for ophthalmic surgery.

Manufacturing materials :

- Polyethylene : needle fitting ;
- Stainless steel: needle tube ;
- Latex : none.

Description :

Sterile, disposable peribulbar, retrobulbar and peri-retrobulbar anaesthetic needles with short 22° bevel (Atkinson).

Instructions for use :

These anesthesia needles are designed for administration of topical anaesthetics in ophthalmic surgery. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the needle aseptically. Do not resterilise. At the end of surgery, place the needle in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g

Construction : Synthetic origin

Sterilisation :

The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

PERIBULBAR ANESTHESIA NEEDLES



R 1646C

23 G - 0,60 x 22 mm

R 1647C

25 G - 0,50 x 22 mm

R 1648C

27 G - 0,40 x 22 mm

PERI-RETROBULBAR ANESTHESIA NEEDLES



R 5107C

23 G - 0,60 x 32 mm

R 5108C

25 G - 0,50 x 32 mm

R 5109C

27 G - 0,40 x 32 mm

RETROBULBAR ANESTHESIA NEEDLES



R 1636C

23 G - 0,60 x 38 mm

R 1637C

25 G - 0,50 x 38 mm

PRODUCT INFORMATION

IRRIGATING CYSTOTOMES

Tip designed for opening of anterior capsule (capsulorhexis) prior to phakoemulsification with a Luer lock to allow irrigation for each cystotome; ready for use in cataract surgery.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Sterile, disposable formed irrigating cystotomes with a 90° angled tip.

Instructions for use :

Irrigating cystotomes are used to perform capsulorhexis. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the cystotome aseptically and place on the operating table. Do not resterilise. At the end of surgery, place the cystotome in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cystotomes. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g

Construction : Synthetic origin

Sterilisation :

The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

IRRIGATING CYSTOTOMES



R 1610C

25 G - 0,50 x 16 mm

R 1611C

27 G - 0,40 x 16 mm



RHEXIS CATARACT SURGERY

PRODUCT INFORMATION HYDRODISSECTION CANNULAS

The tip is designed to separate the cortex from the capsular bag. The flat end allows gentle insertion and suitable irrigation so as to separate the nucleus and cortex from the capsule prior to phakoemulsification.
Ready for use in cataract surgery.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Sterile, disposable hydrodissection cannulas, 45° angled or curved at 8 mm from the tip.

Instructions for use :

Hydrodissection cannulas are designed to dissect the crystalline lens nucleus.
Check that the individual sterile protective pack is unopened and verify the expiry date.
Open the individual sterile protective pack aseptically. Remove the cannula aseptically and place on the operating table.
Do not resterilise. At the end of surgery, place the cannula in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards.
There are two types of waste material : polyethylene and the metal used for the ends of the cannula.
Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g - **Construction :** Synthetic origin

Sterilisation : The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

ANGLED HYDRODISSECTION CANNULA




R 5037C

25 G - 0,50 x 22 mm

R 5039C

27 G - 0,40 x 22 mm

CURVED HYDRODISSECTION CANNULA



R 5038C

27 G - 0,40 x 22 mm



PRODUCT INFORMATION RYCROFT IRRIGATION CANNULAS

Cannulas with very fine angled tips to allow irrigation and holding of the anterior chamber during surgery. The irrigation cannula can be used to provide regular irrigation of the cornea.
Ready for use in cataract surgery.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Sterile, disposable Rycroft anterior chamber irrigation cannula, 40° angled at 4 mm from the tip.

Instructions for use :

Rycroft cannulas are used to irrigate the anterior chamber and moisten the cornea during cataract surgery.
Check that the individual sterile protective pack is unopened and verify the expiry date.
Open the individual sterile protective pack aseptically. Remove the cannula aseptically and place on the operating table.
Do not resterilise. At the end of surgery, place the cannula in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions : Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated : The waste generated is to be considered as waste involving infectious or related hazards.

There are two types of waste material : polyethylene and the metal used for the ends of the cannulas.
Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g - **Construction :** Synthetic origin

Sterilisation : The products are sterilised by irradiation in accordance with standard EN 552.

Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

RYCROFT IRRIGATION CANNULAS



R 1271C

25 G - 0,50 x 22 mm

R 1272C

27 G - 0,40 x 22 mm

R 1273C

30 G - 0,30 x 22 mm



PRODUCT INFORMATION

VISCO-ELASTIC CANNULA

Ready for use for intraocular injection of the viscoelastic solutions.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Disposable, sterile cannula for viscoelastic solutions, 45° angled at 9mm from the tip.

Instructions for use :

Viscoelastic cannulas are used to inject viscoelastic solutions, during eye surgery. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the cannulas aseptically and place on the operating table. Do not resterilise. At the end of surgery, place the cannulas in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g

Construction : Synthetic origin

Sterilisation :

The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

VISCO-ELASTIC CANNULA



R 5011C

25 G - 0,50 x 22 mm



PRODUCT INFORMATION

CORTEX EXTRACTOR (CHARLEUX)

This cannula has a rounded, foam-covered tip and is specially designed for aspirating cortex. Ready to use for cataract surgery.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Sterile, disposable Cortex extractor, curved, with rounded edges tip.

Instructions for use :

Cortex extractor are used to aspirate cortex during cataract surgery. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the cannula aseptically and place on the operating table. Do not resterilise. At the end of surgery, place the cannula in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g - Construction : Synthetic origin

Sterilisation :

The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

CORTEX EXTRACTOR (CHARLEUX)



R 5020C

23 G - 0,60 x 22 mm



PRODUCT INFORMATION

J-SHAPED ASPIRATION CANNULA

The tip of this cannula is specially designed to remove cortex located between 10 o'clock and 2 o'clock. Ready for use in cataract surgery.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Sterile, disposable J-Shaped aspiration cannula.

Instructions for use :

J-Shaped aspiration cannulas are used to aspirate cortex around 12 o'clock. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the cannula aseptically and place on the operating table. Do not resterilise. At the end of surgery, place the cannula in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g - **Construction :** Synthetic origin

Sterilisation : The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

J-SHAPED ASPIRATION CANNULA



R 5003C

25 G - 0,50 x 22 mm



PRODUCT INFORMATION

SILICON TIP POLISHER

Cannula designed for polishing and cleaning the posterior capsule including the equatorial area, prior to insertion of the intraocular lens. Ready for use in cataract surgery.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Silicon : soft polishing sleeve ;
- Latex : none.

Description :

Sterile, disposable Silicon Tip polisher, tube 45° angled at 9 mm from tip, 12 mm soft polishing sleeve.

Instructions for use :

Silicon Tip polishers are used to polish the posterior capsule during cataract surgery. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the Silicon Tip polisher aseptically and place on the operating table. Do not resterilise. At the end of surgery, place the Silicon Tip polisher in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g - **Construction :** Synthetic origin

Sterilisation : the products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

SILICON TIP POLISHER



R 1405C

27 G - 0,40 x 22 mm



PRODUCT INFORMATION

STRAIGHT LACRIMAL CANNULA

Ready to use for cleansing of tear ducts.

Manufacturing materials :

- Polyethylene : cannula fitting ;
- Stainless steel : cannula tube ;
- Latex : none.

Description :

Sterile, disposable straight lacrimal cannula.

Instructions for use :

Lacrimal cannulas are used for cleansing of tear ducts. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the cannula aseptically and place on the operating table. Do not resterilise. At the end of surgery, place the cannula in a sharps disposal box.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 1,80 g

Construction : Synthetic origin

Sterilisation :

the products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

STRAIGHT LACRIMAL CANNULA



R 1205C

23 G - 0,60 x 10 mm



PRODUCT INFORMATION

EYE SHIELD

Transparent anatomically-shaped post-surgery eye protector. Ready for use for cataract surgery.

Manufacturing material :

- Transparent acrylic material ;
- Latex : none.

Description :

Sterile, disposable is.

Instructions for use :

Eye Shield are used to protect the eye after surgery. Check that the individual sterile protective pack is unopened and verify the expiry date. Open the individual sterile protective pack aseptically. Remove the ocular protector aseptically and place on the operating table. Do not resterilise. At the end of surgery, discard in a disposable waste bag.

Packaging :

Individually packaged in sterile blister packs.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

The waste generated is to be considered as waste involving infectious or related hazards. There are two types of waste material : polyethylene and the metal used for the ends of the cannulas. Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 5,23 g

Construction : the different components are of synthetic origin.

Sterilisation :

The products are sterilised by irradiation in accordance with standard EN 552.


Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

EYE SHIELD



R 1040C

Universal transparent ocular shell



PRODUCT INFORMATION

EYES-SPONGES

Designed to ensure high absorption of fluids and blood from the ocular surgical field during cataract surgery.
Ready for use in cataract surgery.

Manufacturing materials :

- PVA (fibre-free, surgical grade) ;
- Polyethylene : stem ;
- Latex : none.

Description :

Sterile, disposable eyes-sponges.

Instructions for use :

The eyes-sponges are designed to absorb fluids at the site of ophthalmic surgery.
Check that the individual sterile protective pack is unopened and verify the expiry date.
Open the individual sterile protective pack aseptically.
Remove the eyes-sponges aseptically and place on the operating table.

The eyes-sponges must be moistened prior to use to get them soft and absorbent. Do not use at the dry state for refractive surgery (Lasik, etc.)

Do not resterilise. At the end of surgery, place the sponges in disposable waste bags.

Packaging :

Individual packs containing 5 sterile micro-sponges.

Storage conditions : Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated : The waste generated is to be considered as waste involving infectious or related hazards.

There are two types of waste material : polyvinylalcohol and polyethylene for the shaft.
Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 0,44 g - **Construction :** the different components are of synthetic origin.

Sterilisation : the products are sterilised by irradiation in accordance with standard EN 552.

Shelf life : the shelf life years is three years.

Manufacturing site :


LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

EYES-SPONGES



R 1061C



PRODUCT INFORMATION

ABSORBENT STICKS

Ready for use in cataract surgery.

Manufacturing materials :

- Cellulose (fibre-free, surgical grade) ;
- Latex : none.

Description :

Sterile, disposable absorbent sticks measuring 5 x 66 mm with bevelled tips.

Instructions for use :

The absorbent sticks used to absorb fluids at the surgical site.

Check that the individual sterile protective pack is unopened and verify the expiry date.

Open the individual sterile protective pack aseptically. Remove all the absorbance sticks from the pack aseptically and place on the operating table.

Do not resterilise. At the end of surgery, place the sticks in disposable waste bags.

Packaging :

Individual packs containing 8 sterile sticks.

Storage conditions :

Store at ambient temperature and protect from moisture.

Disposal conditions :

Waste generated :

Waste generated : The waste generated is to be considered as waste involving infectious or related hazards.

There are two types of waste material : cellulose and polyethylene for the stem.

Disposal conditions comprise incineration taking steps to ensure energy recuperation and environmental protection.

Product weight : 0,29 g.

Construction : the different components are of synthetic origin.

Sterilisation :

The products are sterilised by irradiation in accordance with standard EN 552.

Shelf life : the shelf life years is three years.

Manufacturing site :

LCA SA 28000 CHARTRES France.

CE marking : CE 0120.

ABSORBENT STICKS











R 1060C



RHEXIS CATARACT SURGERY

PRODUCT INFORMATION

SYMBOLS

- CE marking 
- Sterile, by radiations 
- Single use 
- Scratch date 
- Number of batch 
- Which can be recycled 
- Conditioning 
- Unit of order 



REFERENCE	DESIGNATION	PAGE
-----------	-------------	------

ANESTHESIA NEEDLES

R 1636C	Retrobulbar anesthesia needles 38 mm - 23 G - CE 0120	4
R 1637C	Retrobulbar anesthesia needles 38 mm - 25 G - CE 0120	4
R 1646C	Peribulbar anesthesia needles 22 mm - 23 G - CE 0120	4
R 1647C	Peribulbar anesthesia needles 22 mm - 25 G - CE 0120	4
R 1648C	Peribulbar anesthesia needles 22 mm - 27 G - CE 0120	4
R 5107C	Peri-retrobulbar anesthesia needles 32 mm - 23 G - CE 0120	4
R 5108C	Peri-retrobulbar anesthesia needles 32 mm - 25 G - CE 0120	4
R 5109C	Peri-retrobulbar anesthesia needles 32 mm - 27 G - CE 0120	4

CYSTOTOMES

R 1610C	Irrigating cystotomes 25 G - CE 0120	5
R 1611C	Irrigating cystotomes 27 G - CE 0120	5

HYDRODISSECTION CANNULAS

R 5037C	Angled hydrodissection cannula 25 G - CE 0120	6
R 5038C	Curved hydrodissection cannula 27 G - CE 0120	6
R 5039C	Angled hydrodissection cannula 27 G - CE 0120	6

IRRIGATION CANNULAS

R 1271C	Rycroft irrigation cannulas 25 G - CE 0120	7
R 1272C	Rycroft irrigation cannulas 27 G - CE 0120	7
R 1273C	Rycroft irrigation cannulas 30 G - CE 0120	7

VISCO-ELASTIC CANNULA

R 5011C	Visco-elastic cannula 25 G - CE 0120	8
---------	--------------------------------------	---

CORTEX ASPIRATION

R 5020C	Cortex extractor (Charleux) 23 G - CE 0120	9
R 5003C	J-Shaped aspiration cannula 25 G - CE 0120	10
R 1405C	Capsule polisher 27 G - CE 0120	11

LACRIMAL CANNULA

R 1205C	Straight lacrimal cannula 23 G, 10 mm - CE 0120	12
---------	---	----

OCULAR SURGERY ACCESSORY PRODUCTS

R 1040C	Eye shield - CE 0120	13
R 1061C	Eyes-sponges (sachet of 5) - CE 0120	14
R 1060C	Absorbent sticks (sachet of 8) - CE 0120	15