



The future of vision

ALOS EYELID IMPLANTS



ALOS Eyelid Implants

| | |
|-----------------------------|---|
| Brand name | ALOS |
| Models | GOLD (TH) and PLATINUM (PTH) |
| Description and Indications | <ul style="list-style-type: none"> • Biocompatible materials • Due to their fine, spherical design, the implants adapt to the curvature of the eyeball to achieve a better aesthetic result. • Rounded edges to avoid the extrusion of the implant in the long-term. • Suture holes with fixation channels. • Purity of 99,9%(Platinum) and 99,9%(Gold). • Sterile. <p><u>Indication:</u> Eyelid implants are designed for the correction of lagophthalmos.</p> |

GOLD*

| References | TH 06 | TH 08 | TH 10 | TH 12 | TH 14 | TH 16 | TH 18 |
|----------------------|---------|-------|-------|-------|-------|-------|-------|
| Weight (g) | 0.6 | 0.8 | 1.0 | 1.2 | 1.4 | 1.6 | 1.8 |
| Length (mm) | 10.5 | 12.5 | 15.5 | 17.5 | 19.5 | 21.5 | 24.5 |
| Width (mm) | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 |
| Thickness (mm) | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 |
| Sterilization method | Steam. | | | | | | |
| Shelf-life | 3 years | | | | | | |

PLATINUM*

| References | PTH06 | PTH08 | PTH 10 | PTH 12 | PTH 14 | PTH 16 | PTH 18 |
|----------------------|---------|-------|--------|--------|--------|--------|--------|
| Weight (g) | 0.6 | 0.8 | 1.0 | 1.2 | 1.4 | 1.6 | 1.8 |
| Length (mm) | 10.2 | 12.4 | 15.0 | 17.0 | 18.8 | 21.0 | 24.0 |
| Width (mm) | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 |
| Thickness (mm) | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 |
| Sterilization method | Steam. | | | | | | |
| Shelf-life | 3 years | | | | | | |



| | |
|------------|-----------------------------------|
| Brand name | ALOS SIZER SET |
| Reference | SET MEDIDOR PESAS TH |
| Weight (g) | 0.6 - 0.8 - 1.0 - 1.2 - 1.4 - 1.6 |

* Consult other measures

Product with CE marking

REV: 07B. 07/2021