



The future of vision

SYMBLEPHARON RING



Symblepharon Ring

Brand name	AJL SYMBLEPHARON RING		
References	SF-S-AJL	SF-M-AJL	SF-L-AJL
Size	Small: 22 x 19 x 1.2 mm	Medium: 24 x 21 x 1.2 mm	Large: 25 x 22 x 1.2 mm
Description and Indications	AJL SYMBLEPHARON RING is a medical device designed with a geometry specifically adapted to the eyeball. It's used to avoid total or partial adhesion between the internal face of the eyelid and the eyeball, as well as maintaining the depth of the conjunctival bag.		
Material	High resistance acrylic material suitable for ocular use.		
Sterilization method	Ethylene Oxide (ETO).		
Supply	Sterile, single container.		
Shelf-life	4.5 years		

* Product with CE marking