



The future of vision

OCULAR SHIELD



Ocular Shield

<i>Brand name</i>	AJL OCULAR SHIELD		
<i>References</i>	PROCUL-S-AJL	PROCUL-M-AJL	PROCUL-L-AJL
<i>Size</i>	Small: 22 x 19 x 1.7 mm	Medium: 24 x 21 x 1.7 mm	Large: 25 x 22 x 1.7 mm
<i>Description and Indications</i>	AJL OCULAR SHIELD is a medical device specially designed to protect the eyeball. It's designed to protect the eyeball during eyelid surgery.		
<i>Material</i>	High resistance acrylic material suitable for ocular use.		
<i>Sterilization method</i>	Ethylene oxide (ETO).		
<i>Supply</i>	Sterile, single container. Device which includes a bulb syringe, latex free, that facilitates the positioning and extraction of the ocular protector.		
<i>Shelf-life</i>	4.5 years		

* Product with CE marking